

Green Chemistry as a Practical Framework for Pharmaceutical Quality Control and Chemical Laboratories

 Prabakaran Ramachandran^{1*}, Murugesan Palanivelu², Selvaraj Elumalai³, Jayaram Kamma⁴
¹Granules Pharmaceutical Inc, Virginia, USA

²Carnegie Pharmaceutical LLC, New Jersey, USA

³Granules Pharmaceutical Inc, Virginia, USA

⁴Lupin, New Jersey, USA

 DOI: <https://doi.org/10.36348/sjmps.2026.v12i07.005>

| Received: 28.05.2026 | Accepted: 17.07.2026 | Published: 20.07.2026

 *Corresponding author: Prabakaran Ramachandran
 Granules Pharmaceutical Inc, Virginia, USA

Abstract

The application of green chemistry principles in pharmaceutical quality control (QC) and chemical laboratories offers a systematic pathway to reduce environmental impact while maintaining analytical performance. This communication synthesizes key concepts from green analytical chemistry (GAC), including waste minimization, reagent substitution, energy reduction, and the use of standardized greenness metrics such as the Analytical Eco-Scale (AES), AGREE software, and the Green Analytical Procedure Index (GAPI). Furthermore, the transition toward white analytical chemistry (WAC) is discussed as a holistic approach that balances environmental sustainability (green), analytical quality (red), and economic practicality (blue). Collectively, these strategies enable pharmaceutical and chemical laboratories to operate more safely, efficiently, and sustainably without compromising data integrity.

Keywords: Green analytical chemistry; pharmaceutical QC; sustainability metrics; white analytical chemistry; eco-scale; AGREE; GAPI.

Copyright © 2026 The Author(s): This is an open-access article distributed under the terms of the Creative Commons Attribution **4.0 International License (CC BY-NC 4.0)** which permits unrestricted use, distribution, and reproduction in any medium for non-commercial use provided the original author and source are credited.

1. INTRODUCTION

Traditional analytical methods employed in pharmaceutical quality control and chemical testing laboratories are often resource intensive. Common practices involve the use of large solvent volumes, toxic reagents, multi-step sample preparations, and energy-demanding instrumentation. These procedures generate substantial hazardous waste and pose risks to both laboratory personnel and the environment.

Green chemistry, originally formalized by Anastas and Warner (1998), provides a set of twelve design principles aimed at reducing or eliminating hazardous substances. When translated into the analytical domain, these principles form the foundation of green analytical chemistry (GAC). Unlike synthetic chemistry, where the product is a physical substance, analytical chemistry produces information. Therefore, GAC metrics must account for sample handling, energy use, reagent toxicity, and waste generation per analysis.

This communication outlines how pharmaceutical QC and chemical laboratories can adopt green chemistry principles in practical terms, supported by validated assessment tools.

2. Core Benefits of Green Chemistry in Laboratory Settings

2.1 Waste Reduction

Conventional analytical methods generate waste volumes often exceeding 50 mL or 50 g per sample. Green analytical methods aim to reduce this to <1 mL or <1 g per sample. This directly lowers disposal costs and environmental liability.

2.2 Safer Reagents

Reagents classified under the Globally Harmonized System (GHS) with hazard pictograms for toxicity, corrosivity, or flammability are assigned penalty points in metrics such as the Analytical Eco-Scale. Green methods prioritize reagents with NFPA health ratings of 0 or 1 (minimal hazard).

2.3 Miniaturization and Automation

The use of micro-analysis (1–10 mg or μL) and automated, miniaturized instrumentation reduces solvent consumption, analyst exposure, and procedural error. This aligns with GAC principles 2 (minimize sample size) and 5 (automation and miniaturization).

2.4 Energy Efficiency

Instrument energy consumption is categorized as:

- Low: ≤ 0.1 kWh per sample (e.g., FTIR, UV-Vis, UPLC)
- Medium: ≤ 1.5 kWh per sample (e.g., AAS, GC, ICP-MS)
- High: > 1.5 kWh per sample (e.g., NMR, LC-MS, XRD)

Green methods favor low-energy instrumentation where analytically feasible.

2.5 Direct and In-Situ Analysis

Off-line analysis requiring sample collection, preservation, transport, and storage is the least green. In-line or on-line measurements eliminate these steps, reducing energy and material use while improving turnaround time.

3. Tools for Assessing Greenness

Pharmaceutical QC laboratories can objectively evaluate method sustainability using validated metrics. Table 1 summarizes the most widely applied tools.

Table 1: Common green analytical chemistry metrics

Metric	Format	Scoring Principle	Interpretation
Analytical Eco-Scale (AES)	Numerical (0–100)	Start at 100; subtract penalty points for hazards, waste, energy	> 75 = excellent green analysis
AGREE	Circular pictogram + score (0–1)	Algorithm based on 12 GAC principles	Closer to 1.00 = greener
GAPI	Five Pentagrams (15 fields)	Color-coded: green (low impact), yellow (medium), red (high)	Visual overview from sampling to disposal
NEMI	Circular quadrants	Dichotomous (green/white) based on PBTs, hazards, corrosivity, waste volume	Simple but less nuanced

4. Transition to White Analytical Chemistry (WAC)

The limitation of green-only metrics is the exclusion of analytical quality and economic feasibility. White analytical chemistry addresses this by integrating three dimensions:

- Red (Analytical Performance): Accuracy, precision, limit of detection (LOD), linearity.
- Green (Eco-friendliness): Waste, toxicity, energy, renewability.
- Blue (Practical/Economic): Cost per sample, throughput, simplicity, equipment availability.

A method achieving balance across all three dimensions is considered "white" representing an optimal compromise. Tools such as the RGB12 algorithm and the Hexagon metric have been developed to quantify WAC.

5. Practical Recommendations for Pharmaceutical QC Laboratories

Based on the reviewed principles and metrics, the following actions are recommended:

1. Replace toxic solvents (e.g., acetonitrile, methanol) with greener alternatives (e.g., ethanol, water, or supercritical CO_2) where selectivity permits.
2. Reduce sample size to micro or ultra-micro scales (< 10 mg or μL) when consistent with detection limits.
3. Prefer direct analytical techniques (e.g., FTIR, Raman) over methods requiring derivatization.

4. Automate liquid handling and sample preparation to minimize manual steps and reagent volumes.
5. Select low-energy instruments for routine QC applications when high-resolution techniques are not required.
6. Use greenness assessment tools (AES, AGREE, GAPI) during method development and validation.
7. Document waste treatment (recycling, degradation, or passivation) rather than defaulting to "no treatment."

6. CONCLUSIONS

Green chemistry offers a practical, actionable framework for improving the environmental footprint of pharmaceutical QC and chemical laboratories. The transition from traditional to green analytical methods does not require compromise on analytical quality. On the contrary, the adoption of miniaturization, automation, and direct analysis often improves precision, reduces analyst error, and lowers operating costs.

The emergence of white analytical chemistry further refines this approach by ensuring that greenness is not achieved at the expense of functionality or affordability. Routine use of metrics such as AES, AGREE, and GAPI during method development enables objective comparison and continuous improvement.

As regulatory bodies increasingly align with the United Nations Sustainable Development Goals

(particularly SDG 12: Responsible Consumption and Production), the adoption of green and white analytical chemistry is transitioning from a voluntary initiative to a competitive and regulatory necessity.

REFERENCES

- Elga. Type I Ultrapure Water is crucial for HPLC and UHPLC. Available from URL: <https://www.elgalabwater.com/sites/default/files/2018-10/elga-application-note-hplc-and-uhplc%20%281%29.pdf>
- My Green Lab. Autoclaves and Sterilizers. Available from URL: <https://www.mygreenlab.org/autoclaves.html>
- The Royal Society of Chemistry Available from URL: <https://www.rsc.org/>
- The University of Birmingham. Available from URL: <https://www.birmingham.ac.uk/index.aspx>
- Based on the chapter by Palencia M. *et al.*, (2025). Green Analytical Chemistry Metrics. In: Green Environmental Chemical Analysis. Royal Society of Chemistry. DOI: 10.1039/9781837674442-00001
- Merck. Merck Life Science Products & Service Solutions. Available from URL: https://www.sigmaaldrich.com/MT/en?gclid=EAIaIQobChMIicbEsuWWgwMVSJ5oCR2uwQrsEAA YASAAEgKEKvD_B wE
- NFPA. Hazardous materials identification. Available from URL: <https://www.nfpa.org/News-and-Research/Publications-and-media/Blogs-Landing-Page/NFPA-Today/Blog-Posts/2021/11/05/Hazardous-Materials-Identification>
- Repsol. What is energy consumption and why is it important? Available from: <https://www.repsol.com/en/energy-and-the-future/future-of-the-world/energyconsumption/index.cshhtml>
- Green aspects in quality control <https://www.um.edu.mt/library/oar/handle/123456789/123316>
- Greener and Whiter Analytical Chemistry Using Cyrene as a More Sustainable and Eco-Friendlier Mobile Phase Solvent in Chromatography